

Features of the Accu-Chek® Blood Glucose Meters covered by your plan



**Accu-Chek®
Guide Me Meter**



**Accu-Chek®
Guide Meter**

Language	English with universal icons	36 options
Memory	720 test results	720 test results
Battery	Two 3.0-volt lithium coin cell (CR2032) batteries	Two 3.0-volt lithium coin cell (CR2032) batteries
Meter size (inches H x W x D)	3.0 X 1.9 X 0.6	3.2 X 1.9 x 0.8
Test time (seconds)	Less than 4 seconds	Less than 4 seconds
Testing illumination	LCD	LCD with backlight and test strip port light
Meal flagging	On health care apps and smartphones	On meter and health care apps
Minimum blood sample size (µL)	0.6 µL with full-end yellow dosing area	0.6 µL with full-end yellow dosing area
Wireless connectivity	Yes	Yes
Reviewing results	Pre- and post-meal 7, 14, 30 and 90 day averages available on meter or on mySugr™ app	Pre- and post-meal 7, 14, 30 and 90 day averages available on meter or on mySugr™ app
Strip compatibility	Accu-Chek® Guide Test Strips	Accu-Chek® Guide Test Strips

Covered test strips include Accu-Chek® Guide, Accu-Chek® Aviva Plus and Accu-Chek® SmartView. This does not represent complete meter specification information. Visit accu-chek.com/us for complete product specifications and product demonstration videos. Accu-Chek, Accu-Chek Guide Me, Accu-Chek Guide, Accu-Chek Aviva Plus and Accu-Chek SmartView are trademarks of Roche.

Features of the OneTouch® Blood Glucose Meters covered by your plan



OneTouch Ultra® 2 Meter



OneTouch Verio Flex® Meter

Language	English and Spanish	Universal
Memory	500 blood glucose and control solution tests	500 blood glucose and control solution tests
Battery	CR2032 coin cells (two)	CR2032 coin cell
Meter size (inches H x W x D)	3.12 x 2.25 x 0.90	3.38 x 2.05 x 0.63
Test time (seconds)	5	5
Testing illumination	Backlit color display	No
Meal tagging	Add meal tags and comments on meter	Only with the OneTouch Reveal® mobile app
Minimum blood sample size (µL)	1.0 µL	0.4 µL
Reviewing results	<ul style="list-style-type: none"> • See if your results are in or out of range, pre- and post-meal and 7, 14 and 30 day averages • See how food affects your blood sugar • View results on your meter or download the OneTouch Reveal® mobile app 	<ul style="list-style-type: none"> • See if your results are in or out of range • View results on your meter or download the OneTouch Reveal® mobile app
Strip compatibility	OneTouch Ultra® Test Strips	OneTouch Verio® Test Strips

Covered test strips include OneTouch Verio® and OneTouch Ultra.®

This does not represent complete meter specification information. Visit onetouch.com for complete product specifications and product guides. ©2024 LifeScan IP Holdings, LLC. All rights reserved. US-OTB-2200100

Plans are insured through UnitedHealthcare Insurance Company or one of its affiliated companies, a Medicare Advantage organization with a Medicare contract. Enrollment in the plan depends on the plan’s contract renewal with Medicare. Benefits, features and/or devices vary by plan/area. Limitations, exclusions and/or network restrictions may apply.

© 2024 United HealthCare Services, Inc. All Rights Reserved.

Y0066_SPRJ83990_062824_M

SPRJ83990